NextFlow Introduction

Maxime Borry, Quentin Letourneur 07/09/2018

Introduction

The goal of this exercise is to recreate the pipeline showned in Figure 1 represented by its DAG. First clone the GitHub repository and cd into it. Name your pipeline coverage.nf and save it at the root of the cloned repository. You can then execute it with the command nextflow run coverage.nf -with-dag workflow.nf

Channel factory

You may want to use the have a look at the fromFilePairs methods to create the channel when working with paired end data

Mapping

You can use the conda directive to add bowtie2 to your pipeline: conda "bioconda::bowtie2"

The mapping command to use is the following:

bowtie2 -q -1 reads_1.fastq -2 reads_2.fastq -x index_prefix -S output.sam -p nb_cpus

Samtools view

You can use the conda directive to add samtools to your pipeline: conda "bioconda::samtools"

The goal is to convert the sam file into a compressed binary bam file using samtools view

The samtools command to use is the following:

samtools view -S -@ nb_cpus -b -o output.bam input.sam

Samtools sort

You can use the conda directive to add samtools to your pipeline: conda "bioconda::samtools"

The goal is to sort the reads mapped on the reference genome by position

The samtools command to use us the following:

samtools sort -@ nb_cpus -o sorted_output.bam input.bam




Figure 1: Directed Acyclic Graph (DAG) of the desired pipeline. Ellipses represent processes, arrows represents channels

Bedtools

You can use the conda directive to add bedtools to your pipeline: conda "bioconda::bedtools"

The goal is to compute the a position specific coverage of the reference genome provided with the aligned reads

The bedtools command to use is the following:

bedtools genomecov -ibam sorted_input.bam -d > output.gcbout

Coverage Stats

You can use the conda directive to add Python 3.6 and Numpy to your pipeline: conda "python=3.6 numpy" The goal is to compute the coverage (mean and median) statistics from the position specific coverage file. The script to use is located in the bin directory, and therefore is detected by Nextflow automatically.

The command to run this script is the following:

bed2coverage input.gcbout